



D-Glucose (K-GLUHKR / K-GLUHKL) Procedure for ChemWell® 2910 Auto-Analyser

Requirements:

- D-Glucose Assay Kit (K-GLUHKR / K-GLUHKL) (provides ~ 1100 assays / ~ 2200 assays).
- K-GLUHK ChemWell® 2910 assay file.
- Use in association with the D-Glucose Assay Kit (K-GLUHKR / K-GLUHKL) product data booklet.


Use:

For the specific measurement of D-fructose and/or D-glucose especially in wines, fruit juices, beverages and food products.

For specific sample preparation methods refer to the D-Glucose Assay Kit (K-GLUHKR / K-GLUHKL) data booklet.

Assay Principle:

Conversion of D-glucose via the following reactions is directly proportional to the coupled formation of NADPH:



Procedure:

Prepare the assay reagents and calibrators and use with the K-GLUHK ChemWell® 2910 assay file.

D-Glucose Assay Kit Components:

- Bottle 1:** Buffer (25 mL, pH 7.6) plus sodium azide (0.02% w/v) as a preservative. Stable for > 2 years at 4°C.
- Bottle 2:** NADP⁺, ATP and PVP. Stable for > 5 years at -20°C.
- Bottle 3:** Hexokinase plus glucose-6-phosphate dehydrogenase suspension (2.25 mL). Stable for > 2 years at 4°C.

Preparation of Kit Components:

1. Use the contents of bottle 1 as supplied. Stable for > 2 years at 4°C.
2. Dissolve the contents of bottle 2 in 12 mL of distilled water. **Stable for > 1 year at 4°C** or > 2 years at -20°C (to avoid repetitive freeze / thaw cycles, divide into appropriately sized aliquots and store in polypropylene tubes).
3. Use the contents of bottle 3 as supplied. Swirl the bottle to mix contents before use. Stable for > 2 years at 4°C.

Preparation of Assay Reagents: (per ~ 105 assays)

Reagent 1:

Component	Volume
distilled water	19 mL
bottle 1 (buffer)	1 mL
*bottle 2 (NADP/ATP)	1 mL
Total volume	21 mL

*after adding 12 mL of distilled water

Reagent 1 stability: > 30 days at 4°C

Reagent 2:

Component	Volume
distilled water	1.9 mL
bottle 3 (HK/G6P-DH)	0.2 mL
Total volume	2.1 mL

Reagent 2 stability: > 30 days at 4°C

Calibrators:

- K-GLUHK 1: 0 g/L (use distilled water)
K-GLUHK 2: 1.5 g/L D-glucose
K-GLUHK 3: 3.0 g/L D-glucose
K-GLUHK 4: 6.0 g/L D-glucose

Assay Parameters:

- Assay volumes: Reagent 1: 0.200 mL
Sample: 0.003 mL
Reagent 2: 0.020 mL
- Calibrators: 0, 1.5, 3.0, 6.0 g/L D-glucose
Reaction time: 5 min at 37°C
Wavelength: 340 nm
Assay type: endpoint
Reaction direction: increase
Linearity: up to 6 g/L of D-glucose

