

## AMBERSEP 200 H<sup>+</sup> ION EXCHANGE RESIN

02/24

## **G-AMBH**

Ambersep 200 H is a strong acid, cation exchange resin based on sulphonic acid exchange groups on a polystyrene matrix. Its high degree of crosslinking gives stability to the structure of the resin resulting in higher resistance to chemical oxidation and to degradation by mechanical, thermal or osmotic shocks than other comparable cation exchange resins.

**PROPERTIES:** 

Matrix: Styrene divinylbenzene copolymer

Functional group: Sulphonic acid (SO<sub>3</sub><sup>-</sup>)

Physical form: Beige coloured beads

Ionic form:

**Total exchange capacity:**  $\geq 1.6 \text{ meq/mL (H}^+)$ 

**Moisture capacity:** 50 - 57% (H<sup>+</sup>)

Shipping weight: 770 g/L

 Harmonic mean size:
 0.82 to 1.00 mm

 Fines:
 < 0.60 mm 1.0% max</td>

 Oversize:
 > 1.18 mm 15% max

Reversible swelling (max):  $Na^+ \rightarrow H^+ 6\%$ 

## **SUGGESTED OPERATING CONDITIONS:**

Max. temperature: 135°C

Working flow rate: 10 - 120 Bed Volumes per hour (BV/h)

Regenerant:HCI $H_2SO_4$ Concentration (%):5-61.5-4Regenerant level (g/L):80-200125-250Regeneration flow rate:4-5 (BV/h)4-12 (BV/h)

Minimum contact time: 30 min

Slow rinse: 2 BV at regeneration flow rate Fast rinse: 2 - 4 BV at working flow rate