




WHEAT ARABINOXYLAN (Low Viscosity) (Lot 160419b)

P-WAXYL

05/20

CAS: 9040-27-1

STRUCTURE



Schematic representation of wheat arabinoxylan unit (Ara : Xyl = 38 / 62)

PROPERTIES

Purity:	~ 95% (dw basis)
Sugar Ratio:	Arabinose:Xylose = 38 / 62
Viscosity:	13 cSt (1% w/v; Ostwald C-type viscometer, 30°C)
Molecular Weight:	56.7 Kd (MAALS)
Starch Content:	< 0.05%
Beta-Glucan:	< 0.05%
Protein:	1.0%
Moisture:	4.6%
Ash:	2.0%
Physical Description:	Slightly off-white, odourless powder


STORAGE CONDITIONS

Store dry at room temperature in a well sealed container. Under these conditions, the product is stable for several years.

METHOD OF DISSOLUTION (FOR 1% W/V SOLUTION)

Accurately weigh 1 g of arabinoxylan into a 120 mL dry pyrex beaker. Wet the sample with 8 mL of 95% ethanol. Add a magnetic stirrer bar, followed by 90 mL of distilled water. Immediately place the beaker containing the slurry on a magnetic stirrer-hotplate and heat at a setting of 100°C with vigorous stirring. Loosely cover the beaker with aluminium foil and stir and boil the contents until the arabinoxylan completely dissolves (approx. 10 min). Allow the solution to cool to room temperature with continued stirring. Adjust the volume to 100 mL. The solution may be very slightly opalescent due to the presence of trace amounts of protein. Arabinoxylan solutions can be stored at room temperature for several weeks in a well sealed storage bottle. Microbial contamination is prevented by adding a few drops of toluene to the storage bottle.

Gas liquid chromatography of the alditol acetates derived from hydrolysis and derivatisation of low viscosity wheat arabinoxylan (Lot I60419b)



GLC

A typical polysaccharide sample (~ 10 mg) was hydrolysed using 2N TFA at 120°C for 60 min. Subsequent sodium borohydride reduction was performed in 1N NH₄OH for 90 minutes at 40°C. The corresponding alditol acetates were prepared using acetic anhydride and 1-methyl imidazole, extracted into DCM and analysed by GC. Chromatography was performed on a Shimadzu GC-14B with CHROMATOPACK C-R8A using a Packed glass column (6 ft x 5 mm OD, 3 mm ID) with 3% Silar 10C on W-HP (80-100 mesh). The carrier gas was nitrogen at 130 KPa. Injector temperature; 250°C; Column temperature; 230°C. Detection by FID with 60KPa H₂ pressure and 50 KPa air pressure.